

ADAPTERS

15

# ORAL ROUTE

**DOS'UP®** adapters are used with **DOS'UP®** pipettes and dosing systems. They serve to optimise dosing with improved hygiene.

#### PRODUCT INFORMATION

### Adapters

Material	PE (Food and/or pharmacopoeia grade)
Use	For PP28 bottleneck, glass and plastic
	> Models for 1-hole pipette with Rounded or Connect tip

## STANDARD REFERENCE(S)

ADAPTER FOR PIPETTE					
ltem code	Bottleneck diameter	Bottle	Shape	Pipette	
1DAD0202	PP28	Glass	Rounded	1 hole	
1DAD0203	PP28	Plastic	Rounded	1 hole	
20034500*	PP28	Glass	Connect	1 hole	

### **FOR 1-HOLE PIPETTE**



5AA adapter
For pipette with rounded tip
and PP28 bottleneck



5CA adapter
For connect pipette
and PP28 bottleneck

### **Union Plastic**

ZA Robert BP 27 43140 saint-Didier-en-Velay-France www.union-plastic.com Tel.: +33 **(0)4 71 61 13 09** union-plastic@omerin.com



The information given in this technical data sheet is indicative and subject to change without notice. Union Plastic can under no circumstances be held responsible for any incidents arising from inappropriate use, failure to comply with the rules of the trade and standards in force, or failure to comply with the instructions for use supplied. Can be used without special training or skills. Non-contractual illustrations and photos. Reproduction prohibited without prior consent from Union Plastic.