PM-CARE® REMOV-KIT

IUD REMOVAL KIT









The PM-CARE® REMOV-KIT IUD removal kit is used to remove intrauterine devices from the uterine cavity. It is composed of a Pozzi forceps, a removal forceps, a pair of string scissors and a hysterometer.

IUD REMOVAL KIT COMPONENTS:

- The straight Pozzi forceps grasps the cervix, performing "traction-fixation" on the cervix to enable the flexion angle of the uterus to be adjusted to 180°. This correction of the flexion angle is necessary in order to remove the intrauterine device without risking uterine perforation.
- The IUD removal forceps gently remove and cleanse the uterus and working area.
- Distal string scissors are used to cut the strings of the intrauterine device.
- The hysterometer is used to assess the depth of the uterine cavity.

AVAILABLE REFERENCE(S)

Reference	POZZI FORCEPS Model	REMOVAL FORCEPS Furnished	STRING SCISSORS Model	HYSTEROMETER Diameter (Fr)
PS3620000	Straight	Yes	Distal	14

ABOUT THIS MEDICAL DEVICE

2011
Yes
ls
Yes
Ethylene oxide sterilisation
No
5 years
CE
Box of 10 units in individual sterile wrapping
In its original box in a cool and dry place

PRINCE MEDICAL SAS ZA La Sente Du Moulin 64 rue Gutenberg - 60530 Ercuis - France www.prince-medical.com Tel: +33 (0)3 44 26 54 76 prince-medical@omerin.com



The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner