

# PM-CARE® REMOV

14

## IUD REMOVAL FORCEPS



The PM-CARE® REMOV IUD removal forceps are used during a gynaecological examination

- When removing an intrauterine device. It holds the strings of the device to pull it down for removal.
- When inserting a new device. It makes it possible to hold a compress for disinfecting and cleaning the work area and/or the cervix soiled by the removal and to prepare the insertion of the new device.

### AVAILABLE REFERENCE(S)

Reference	Length (cm)
PS3570025	25

### ABOUT THIS MEDICAL DEVICE

1 <sup>st</sup> release date	2011
Single use	Yes
Class	Is
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 20 units in individual sterile wrapping
Storage	In its original box in a cool and dry place



PRINCE MEDICAL SAS  
 ZA La Sente Du Moulin  
 64 rue Gutenberg - 60530 Ercuis - France  
[www.prince-medical.com](http://www.prince-medical.com)  
 Tel: +33 (0)3 44 26 54 76  
[prince-medical@omerin.com](mailto:prince-medical@omerin.com)

**PRINCE MEDICAL**  
 DEVICES

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.