
# **PM-CARE®**RETRIEVAL

# IUD HOOK



The **PM-CARE**<sup>\*</sup> **RETRIEVAL** IUD hook is used to retrieve IUD strings that have been lost in the endocervical canal or the uterine cavity, as well to remove the IUD, if necessary.

## AVAILABLE REFERENCE(S)



### ABOUT THIS MEDICAL DEVICE

1 <sup>st</sup> release date	2014
Single use	Yes
Class	ls
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 10 units in individual sterile wrapping
Storage	In its original box in a cool and dry place

#### PRINCE MEDICAL SAS

ZA La Sente Du Moulin 64 rue Gutenberg - 60530 Ercuis - France www.prince-medical.com Tel: +33 **(0)3 44 26 54 76** prince-medical@omerin.com



The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

\* Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.