

# PM-CARE® SPECULUM

24

CUSCO OR COLLIN SPECULUM



The **PM-CARE® SPECULUM** speculum is used to explore, view, and examine the uterine cervix and vagina, and to access the cervical canal and uterine cavity when performing various examinations and surgical procedures. The Collin speculum enables larger devices to be inserted. The **PM-CARE® SPECULUM** is made from very high quality medical thermoplastic material.

## AVAILABLE REFERENCE(S)

Reference	Model	Size	Diameter (mm)	Colour
PS3592000	Cusco	XS	20	White
PS3592400	Cusco	S	24	Blue
PS3592600	Cusco	M	26	Red
PS3593000	Cusco	L	30	Green
PS3932700	Collin	S	27	White
PS3933000	Collin	M	33	Green

## ABOUT THIS MEDICAL DEVICE

1 <sup>st</sup> release date	Cusco Speculum: 2012 Collin Speculum: 2015
Single use	Yes
Class	Is
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	• Cusco Speculum: box of 50 units in individual sterile wrapping • Collin Speculum: box of 35 units in individual sterile wrapping
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS  
 ZA La Sente Du Moulin  
 64 rue Gutenberg - 60530 Ercuis - France  
[www.prince-medical.com](http://www.prince-medical.com)  
 Tel: +33 (0)3 44 26 54 76  
[prince-medical@omerin.com](mailto:prince-medical@omerin.com)

**PRINCE MEDICAL**  
 DEVICES

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.