


CLEANBRUSH® KIT-2

60

2-PART CLEANING SWAB KIT



The CLEANBRUSH® KIT-2 swab kits are intended for cleaning the working channels, valves and valve cages of flexible endoscopes.

They are suitable for use with most major brands of gastroscopes, colonoscopes, duodenoscopes, enteroscopes, bronchoscopes, cystoscopes, ureteroscopes, nephroscopes, hysteroscopes and nasolaryngoscopes.

KIT COMPONENTS (SEE THE TABLE FOR THE DETAILED COMPONENTS BY REFERENCE)

- CLEANBRUSH® ENDO: cleaning brush for working channel
- CLEANBRUSH® PISTON: cleaning brush for valves, valve cages
- CLEANBRUSH® EXO-B: brush for external cleaning
- CLEANBRUSH® EXO- S: sponge for external cleaning

ABOUT THIS MEDICAL DEVICE

1 st release date	1999
Single use	Yes
Class	I
Sterile	No (ethylene oxide decontamination)
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 25 or 50 units depending on product reference
Storage	In its original box in a cool and dry place



PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

PRINCE MEDICAL
 DEVICES

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

CLEANBRUSH® KIT-2

61

AVAILABLE REFERENCE(S)

REFERENCES FOR CLEANING BRUSHES

Kit reference	For working channel	For valves	Packaging
AN1850110	AN1050110 (with 1 brush Ø 1.8 mm L100 cm)	AN1061110 (with 2 brushes Ø11 & 5mm, L15.5 cm)	50
AN1853010	AN1053010 (with 1 brush Ø 3 mm L100 cm)		
AN1853018	AN1053018 (with 1 brush Ø 3 mm L180 cm)		
AN1855010	AN1055010 (with 1 brush Ø 5 mm L100 cm)		
AN1855018	AN1055018 (with 1 brush Ø 5 mm L180 cm)		
AN1855020	AN1055020 (with 1 brush Ø 5 mm L200 cm)		
AN1855023	AN1055023 (with 1 brush Ø 5 mm L230 cm)		
AN1856523	AN1056523 (with 1 brush Ø 6.5 mm L230 cm)		
AN1860110	AN1060110 (with 2 brushes Ø 1.8 mm L100 cm)		
AN1863010	AN1063010 (with 2 brushes Ø 3 mm L100 cm)		
AN1863018	AN1063018 (with 2 brushes Ø 3 mm L180 cm)		
AN1865010	AN1065010 (with 2 brushes Ø 5 mm L100 cm)		
AN1865018	AN1065018 (with 2 brushes Ø 5 mm L180 cm)		
AN1865020	AN1065020 (with 2 brushes Ø 5 mm L200 cm)		
AN1865023	AN1065023 (with 2 brushes Ø 5 mm L230 cm)		
AN1866523	AN1066523 (with 2 brushes Ø 6.5 mm L230 cm)		

KIT COMPONENT REFERENCES

Kit reference	Brush for external cleaning	Sponge for external cleaning	Packaging
AN2750000	AN2580000 (U-shaped with soft hair, useful length 2 x 8 cm)	AN1790000	25

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

PRINCE MEDICAL
 DEVICES