


CLEANBRUSH® KIT-3

62

3-PART CLEANING SWAB KIT



The **CLEANBRUSH® KIT-3** swab kits are intended for cleaning the working channels, valves and valve cages of flexible endoscopes.

They are suitable for use with the main gastroscopes, colonoscopes, duodenoscopes, enteroscopes, bronchoscopes, cystoscopes, ureteroscopes, nephroscopes, hysteroscopes and naso-laryngoscopes on the market.

KIT COMPONENTS (SEE THE TABLE FOR THE DETAILED COMPONENTS BY REFERENCE)

- **CLEANBRUSH® ENDO:** cleaning brush for working channel
- **CLEANBRUSH® ENDO (AN1055002, AN1056523):** cleaning brush for aspiration channel
- **CLEANBRUSH® ENDO (AN1050123):** cleaning brush for air/water channel, and/or Water Jet and/or balloon inflation channel
- **CLEANBRUSH® PISTON:** cleaning brush for valves, valve cages
- **CLEANBRUSH® EXO-B:** brush for external cleaning

ABOUT THIS MEDICAL DEVICE

1 st release date	1999
Single use	Yes
Class	I
Sterile	No (ethylene oxide decontamination)
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 25 or 50 units depending on product reference
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

PRINCE MEDICAL
 DEVICES

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

CLEANBRUSH® KIT-3

AVAILABLE REFERENCE(S)

REFERENCES FOR CLEANING BRUSHES

OTHER KIT COMPONENTS

Kit reference	For working channel	For aspiration channel	For valves	For air/water channel, water jet, balloon inflation channel	External cleaning brush	Packaging	
AN1883018	AN1053018 (with 1 brush Ø 3 mm L180 cm)				AN1050123 (with 1 brush Ø 1.8 mm L230 cm)	50	
AN1885023	AN1055023 (with 1 brush Ø 5 mm L230 cm)						
AN1886523	AN1056523 (with 1 brush Ø 6.5 mm L230 cm)						
AN1893018	AN1063018 (with 2 brushes Ø 3 mm L180 cm)						
AN1895023	AN1065023 (with 2 brushes Ø 5 mm L230 cm)						
AN1896523	AN1066523 (with 2 brushes Ø 6.5 mm L230 cm)						
AN2770110	AN1050110 (with 1 brush Ø 1.8 mm L100 cm)						
AN2780110	AN1060110 (with 2 brushes Ø 1.8 mm L100 cm)						
AN2183010	AN1053010 (with 1 brush Ø 3 mm L100 cm)						AN1055002 (with 1 brush Ø 5 mm, L20 cm)
AN2183018	AN1053018 (with 1 brush Ø 3 mm L180 cm)						
AN2193010	AN1063010 (with 2 brushes Ø 3 mm L100 cm)	AN1061110 (with 2 brushes Ø 11 & 5 mm, L15.5 cm)					
AN2193018	AN1063018 (with 2 brushes Ø 3 mm L180 cm)						
AN2953018	AN1053018 (with 1 brush Ø 3 mm L180 cm)		AN1055023 (with 1 brush Ø 5 mm, L230 cm)				
AN2963018	AN1063018 (with 2 brushes Ø 3 mm L180 cm)						
AN4475018	AN1055018 (with 1 brush Ø 5 mm L180 cm)						
AN4475023	AN1055023 (with 1 brush Ø 5 mm L230 cm)						
AN4495018	AN1065018 (with 2 brushes Ø 5 mm L180 cm)		AN1056523 (with 1 brush Ø 6.5 mm, L230 cm)				
AN4495023	AN1065023 (with 2 brushes Ø 5 mm L230 cm)						
AN2520110	AN1050110 (with 1 brush Ø 1.8 mm L100 cm)						
AN2523007	AN1053007 (with 1 brush Ø 3 mm L75 cm)						
AN2523010	AN1053010 (with 1 brush Ø 3 mm L100 cm)						
AN2523018	AN1053018 (with 1 brush Ø 3 mm L180 cm)						
AN2525010	AN1055010 (with 1 brush Ø 5 mm L100 cm)						
AN2525023	AN1055023 (with 1 brush Ø 5 mm L230 cm)						
AN2525030	AN1055030 (with 1 brush Ø 5 mm L300 cm)						
AN2530110	AN1060110 (with 2 brushes Ø 1.8 mm L100 cm)						
AN2533010	AN1063010 (with 2 brushes Ø 3 mm L100 cm)						
AN2533018	AN1063018 (with 2 brushes Ø 3 mm L180 cm)						
AN2535010	AN1065010 (with 2 brushes Ø 5 mm L100 cm)						
AN2535023	AN1065023 (with 2 brushes Ø 5 mm L230 cm)						
AN2535030	AN1065030 (with 2 brushes Ø 5 mm L300 cm)						
						AN2580000 (U-shaped with soft hair, useful length 2 x 8 cm)	25

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.