

CLEANBRUSH® TRACH

66


BRUSH FOR TRACHEAL CANNULA



The CLEANBRUSH® TRACH tracheal cannula brushes are intended for cleaning tracheostomy cannulas.

AVAILABLE REFERENCE(S)

Reference	Brush diameter (mm)	Sheath diameter (mm)	Length (cm)	Packaging
AS2201002	10	3.3	23	25
AS2201402	14	3.3	23	25
AS2201802	18	3.3	23	20



ABOUT THIS MEDICAL DEVICE

1 st release date	2007
Single use	Yes
Class	Is
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 20 or 25 units in sterile wrapping depending on product reference
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

PRINCE MEDICAL
 DEVICES

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.