

ECHO HOOD® APERTURE-1

44

BALLOON 1 OPENING



The ECHO HOOD® APERTURE-1 balloon for echo-endoscope fits on the ultrasound probe to make contact with the gastrointestinal wall and thus provide a good ultrasound image.

ADVANTAGES

- Can withstand handling.
- Adequate thickness for the passage of ultrasound waves.
- Product compatible for use with the main endoscopes on the market.

AVAILABLE REFERENCE(S)

Reference	Compatible endoscopes	Balloon length (mm)	Opening size (mm)
AS1470000	ECHO-DUODENOSCOPE Pentax	34	8.7
AS1830000	ECHO-DUODENOSCOPE Olympus	15	6.4
AS3680000	ECHO-DUODENOSCOPE Olympus 180	15	6.5
AS2800000	ECHO-ENDOSCOPE Fujinon convex	18	5.8
AS3560000	ECHO-ENDOSCOPE Fujinon convex	18	6

ABOUT THIS MEDICAL DEVICE

1 st release date	1998
Single use	Yes
Class	Is
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	Yes
Lifetime	2 years
Certification	CE
Packaging	Box of 20 units in individual sterile wrapping
Storage	In its original box, protected from ultraviolet rays, in a cool and dry place

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

PRINCE MEDICAL
 DEVICES

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.