
# **ENDOLINE® 1V - JECT**



Tel: +33 (0)3 44 26 54 76 prince-medical@omerin.com

26



The ENDOLINE® 1V-JECT 1-way injection needle has been designed specifically and exclusively for haemostasis of gastrointestinal haemorrhages, sclerosis of oesophageal varices and inflating and mark-

The needle is only suitable for use with flexible endoscopes.

- Radiopaque distal part visible on the screen, for models with a metal tip.

### 2.0

Evolution of the handle 1.0, thicker for a better grip in the hand.

## NEW 3.0

**Optimal ergonomics** and lighter weight, this handle has been designed with the user and the environment in mind!

#### ABOUT THIS MEDICAL DEVICE

1 <sup>st</sup> release d	late	1996
Single	use	Yes
CL	ass	lla
Ste	rile	Yes
Sterilisation meth	hod	Ethylene oxide sterilisation
Contains la	itex	No
Lifeti	ime	5 years
Certificat	tion	CE
Packag	ging	Box of 5 units in individual sterile wrapping
Stor	age	In its original box in a cool and dry place

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

### DIGESTIVE & BRONCHIAL ENDOSCOPY/HAEMOSTASIS

# **ENDOLINE® 1V – JECT**

		Reference	Catheter diameter	Needle diameter	<b>Useful length</b>	Needle outlet length
		(mm/Fr)	(mm)	(cm)	(mm)	
HANDLE 1.0	> Metal Tip					
	Piletat rip	AS1990718	2.35/7	0.7	180	5
		AS1990723	2.35/7	0.7	230	5
		AS1991018	2.35/7	1	180	5.5
		AS1991023	2.35/7	1	230	5.5
HANDLE 2.0	> Plastic Tip					
		AS2940518	1.75/5.3	0.5	180	4
		AS2940523	1.75/5.3	0.5	230	4
		AS2940718	2.35/7	0.7	180	5
		AS2940723	2.35/7	0.7	230	5
		AS2941018	2.35/7	1	180	5.5
		AS2941023	2.35/7	1	230	5.5
	> Metal Tip					
		AS3010718	2.35/7	0.7	180	5
		AS3010723	2.35/7	0.7	230	5
		AS3011018	2.35/7	1	180	5.5
		AS3011023	2.35/7	1	230	5.5
		AS3720523	2.35/7	0.5	230	3
		AS4310523	2.35/7	0.5	230	4.5
HANDLE 3.0	> Plastic Tip					
		AS4480518	1.75/5.3	0.5	180	4
		AS4480523	1.75/5.3	0.5	230	4
		AS4480718	2.35/7	0.7	180	5
		AS4480723	2.35/7	0.7	230	5
		AS4481018	2.35/7	1	180	5

#### PRINCE MEDICAL SAS

ZA La Sente Du Moulin 64 rue Gutenberg - 60530 Ercuis - France www.prince-medical.com Tel: +33 **(0)3 44 26 54 76** prince-medical@omerin.com

PRINCE MEDICAL

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

2.35/7

230

1

5

\* Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

AS4481023