


# ENDOLINE® NASO-CONNECT

LUER-LOCK CONNECTOR

21



The ENDOLINE® NASO-CONNECT Luer-Lock connector enables the aspiration of bile or pancreatic juices. It also allows to remove the endoscope while leaving the catheter in place. This component is part of the ENDOLINE® NASO-KIT nasobiliary kit (for more information on the kit, please consult the data sheet page 18).

## AVAILABLE REFERENCE(S)

Reference	Diameter (Fr)	Colour
AS0050500	5	grey
AS0050600	6	blue
AS0050700	7	orange
AS0051000	10	white



## ABOUT THIS MEDICAL DEVICE

1 <sup>st</sup> release date	1997
Single use	Yes
Class	Is
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 50 units in individual sterile wrapping
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS  
 ZA La Sente Du Moulin  
 64 rue Gutenberg - 60530 Ercuis - France  
[www.prince-medical.com](http://www.prince-medical.com)  
 Tel: +33 (0)3 44 26 54 76  
[prince-medical@omerin.com](mailto:prince-medical@omerin.com)

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

PRINCE MEDICAL  
 DEVICES