


# ENDOLINE® NASO-TUBE

NASAL TUBE

23



The ENDOLINE® NASO-TUBE nasal tube is used to retrieve a nasobiliary catheter from the larynx and remove it through the nose. This component is part of the ENDOLINE® NASO-KIT nasobiliary kit (for more information on the kit, please consult the data sheet page 18).

## AVAILABLE REFERENCE(S)

Reference	Nasal tube Int. diam. / Ext. diam. (mm / mm)	Length of the nasal tube (cm)	Compatible with nasobiliary catheter diameter (Fr)
AS0220550	2.0 / 3.0	50	5
AS0220750	2.6 / 3.7	50	7
AS0221050	3.7 / 4.9	50	10



## ABOUT THIS MEDICAL DEVICE

1 <sup>st</sup> release date	1996
Single use	Yes
Class	Is
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 10 units in individual sterile wrapping
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS  
 ZA La Sente Du Moulin  
 64 rue Gutenberg - 60530 Ercuis - France  
[www.prince-medical.com](http://www.prince-medical.com)  
 Tel: +33 (0)3 44 26 54 76  
[prince-medical@omerin.com](mailto:prince-medical@omerin.com)

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

PRINCE MEDICAL  
 DEVICES