


ENDOLINE® SPRAY

CATHETER SPRAY

33



The **ENDOLINE® SPRAY** catheter spray enables staining products to be injected into the digestive system's mucous membrane. This injection is used for mucosal analysis in case of lesions or polyps in the digestive system, in patients with a relevant medical history or for diagnosis purposes.

ADVANTAGES


- Easy one-handed operation.
- Product compatible for use with the main endoscopes on the market.
- Luer-Lock connector for connecting a syringe for injection of staining products.

AVAILABLE REFERENCE(S)

Reference	Catheter diameter (mm)	Useful length (cm)
AS2791823	1.8	230

ABOUT THIS MEDICAL DEVICE

1 st release date	2009
Single use	Yes
Class	Is
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 5 units in individual sterile wrapping
Storage	In its original box in a cool and dry place



PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Erceuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

PRINCE MEDICAL
 DEVICES