

# ENDOLINE® VALVE

50

## WORKING CHANNEL VALVE



VALVE WITH LATERAL IRRIGATION FITTING



VALVE WITH LATERAL IRRIGATION FITTING AND NON-RETURN VALVE

The **ENDOLINE® VALVE** working channel valve is a phosphorescent valve that fits on the endoscope's working channel. It gives accessories access to the working channel, guarantees a tight seal during the examination and eliminates the risk of cross-contamination between patients.

### ADVANTAGES

- Narrow closure to prevent leakage.
- Product compatible for use with the main endoscopes on the market.
- References with trans-valve fitting connecting to commonly used syringe models.
- References with lateral irrigation fitting for lavage through the working channel while maintaining the accessory in place.
- References with integrated non-return valve to prevent gastric reflux in the endoscope.

### ABOUT THIS MEDICAL DEVICE

1 <sup>st</sup> release date	2007
Single use	Yes
Class	Is
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 50 or 100 units in individual sterile wrapping, depending on product reference
Storage	In its original box in a cool and dry place



VALVE ONLY

VALVE WITH FITTING

PRINCE MEDICAL SAS  
 ZA La Sente Du Moulin  
 64 rue Gutenberg - 60530 Ercuis - France  
[www.prince-medical.com](http://www.prince-medical.com)  
 Tel: +33 (0)3 44 26 54 76  
[prince-medical@omerin.com](mailto:prince-medical@omerin.com)



The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

## ECHO HOOD® VALVE

## DIGESTIVE &amp; BRONCHIAL

Reference	Compatible endoscopes	Supplied with a connector	Colour	Packaging
AS2430000	Olympus, Fujinon	-	orange	100
AS2440000	Pentax	-	yellow	100
AS2500000	Olympus, Fujinon	trans-valve	orange	100
AS2510000	Pentax	trans-valve	yellow	100
AS2450000	Olympus, Fujinon	lateral irrigation	orange	50
AS2460000	Pentax	lateral irrigation	yellow	50
AS3810000	Olympus, Fujinon	lateral irrigation & AR valve	orange	50
AS3820000	Pentax	lateral irrigation AR & valve	yellow	50

## BRONCHIAL ENDOSCOPY

AS2690000	Olympus, Fujinon	-	blue	100
AS3020000	Pentax	-	pink	100
AS3350000	Olympus, Fujinon	trans-valve	blue	100
AS3400000	Pentax	trans-valve	pink	100

PRINCE MEDICAL SAS  
 ZA La Sente Du Moulin  
 64 rue Gutenberg - 60530 Ercuis - France  
[www.prince-medical.com](http://www.prince-medical.com)  
 Tel: +33 (0)3 44 26 54 76  
[prince-medical@omerin.com](mailto:prince-medical@omerin.com)

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.