ENDOLINE® Y-TRACT

16









The ENDOLINE® Y-TRACT Y-clamp fitting provides practitioners with 2 injection ports on a balloon catheter with 1 injection port, and also enables the fixation of the guide wire after catheter placement.

AVAILABLE REFERENCE(S)

Reference Description

AS1130300 CH3 Y-clamp fitting

ABOUT THIS MEDICAL DEVICE

1 st release date	1998
Single use	Yes
Class	Is
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 50 units in individual sterile wrapping
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS ZA La Sente Du Moulin 64 rue Gutenberg - 60530 Ercuis - France www.prince-medical.com Tel: +33 (0)3 44 26 54 76 prince-medical@omerin.com



The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner