# PM-LIFE® OVO-KIT

## OOCYTE PUNCTURE KIT









The PM-LIFE® OVO-KIT oocyte puncture kit is used to perform an ultrasound-guided puncture of the ovarian follicles for the purpose of an in vitro fertilisation. The reference PS4341533 also contains a stopcock (with its male Luer-lock connection base) enabling follicular rinsing.

#### **ADVANTAGES**

• The needle is stainless steel with a triple rear bevel and two echogenic areas to capture the ultrasound during ultrasound guidance.

#### **SET COMPONENTS:**

- · A needle and its attachment base (with tube cover).
- 2 or 3 tubes with their Luer Lock connection bases.
- 1 Falcon-type cap.
- · A Stopcock (only for the reference PS4341533).

### AVAILABLE REFERENCE(S)

		NEEDLE		
Reference	Inner diameter (mm)	Outer diameter (mm / Gauge)	Useful length (cm)	Stopcock
PS3251533	1.2	1.5 / 17	31.5	-
PS4341533	1.2	1.5 / 17	31.5	Yes

#### ABOUT THIS MEDICAL DEVICE

1 <sup>st</sup> release date	2010	
Single use	Yes	
Class	lla	
Sterile	Yes	
Sterilisation method	Sterilised via gamma radiation	
Contains latex	No	
Lifetime	5 years	
Certification	CE	
Packaging	PS3251533: box of 5 kits in individual sterile wrapping	
	<ul> <li>PS4341533: box of 4 kits in individual sterile wrapping</li> </ul>	
Storage	In its original box in a cool and dry place	

PRINCE MEDICAL SAS ZA La Sente Du Moulin 64 rue Gutenberg - 60530 Ercuis - France www.prince-medical.com Tel: +33 (0)3 44 26 54 76 prince-medical@omerin.com



The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner

Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.