

OXYCAP® EVO

52

MOUTHPIECE



MOUTHPIECE
WITH 2 CONNECTORS
MOUNTED



MOUTHPIECE
WITH BAND
AND 2 CONNECTORS
MOUNTED



The OXYCAP® EVO mouthpiece is used in upper digestive endoscopies (gastroscopy, duodenoscopy, etc.) to protect the patient's teeth, as well as to protect endoscopes passing through the oral cavity. The connectors allow continuous oxygenation during the procedure.

MOUTHPIECES

Reference	Configuration	Packaging
AS420000	Only	75
AS421000	With band	75
AS425000	With 2 connectors mounted	75
AS426000	With band and 2 connectors mounted	75

ADDITIONAL ACCESSORIES

Reference	Configuration	Packaging
AS422000	Band only	75
AS427000	Connector with Luer	10
AS428000	Connector with filter	10
AS423000	White connector only	100
AS439000	Yellow connector only	100

ABOUT THIS MEDICAL DEVICE

1 st release date	2017
Single use	Yes
Class	Is
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	No
Lifetime	2 years
Certification	CE
Packaging	Box of 10, 75 or 100 units in individual sterile wrapping depending on product reference
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS
ZA La Sente Du Moulin
64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

PRINCE MEDICAL
DEVICES

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.