#### 5

# PM-CARE®INSERT-KIT

# IUD INSERTION KIT







HYSTEROMETER

The **PM-CARE**° **INSERT-KIT** IUD insertion kit is used to insert an intrauterine device in the intrauterine cavity. It is composed of a Pozzi forceps, a Cheron forceps, a pair of string scissors and a hysterometer.

### **IUD INSERTION KIT COMPONENTS:**

- The straight Pozzi forceps grasps the cervix, performing "traction-fixation" on the cervix to enable the flexion angle of the uterus to be adjusted to 180°. This correction of the flexion angle is necessary in order to insert instruments and fit the intrauterine device without risking uterine perforation.
- The angled Cheron forceps cleans the uterus and the working area.
- Lateral or distal string scissors are used to cut the strings of the intrauterine device.

  The practitioner can choose to use a pair of lateral or distal string scissors according to their personal preference and ease in the act to be performed.
- The hysterometer is used to assess the depth of the uterine cavity.

## AVAILABLE REFERENCE(S)

	POZZI FORCEPS	CHERON FORCEPS	STRING SCISSORS	HYSTEROMETER
Reference	Model	Model	Model	Diameter (Fr)
PS4190000	Straight	Angled	Distal	(not supplied)
PS3840000	Straight	Angled	Distal	14
PS3230000	Straight	Angled	Lateral	14

#### ABOUT THIS MEDICAL DEVICE

1st release date	2011	
Single use	Yes	
Class	ls	
Sterile	Yes	
Sterilisation method	Ethylene oxide sterilisation	
Contains latex	No	
Lifetime	5 years	
Certification	CE	
Packaging	Box of 10 units in individual sterile wrapping	
Storage	In its original box in a cool and dry place	

PRINCE MEDICAL SAS
ZA La Sente Du Moulin
64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.



