


CLEANBRUSH® KIT-2

60

2-PART CLEANING SWAB KIT



The **CLEANBRUSH® KIT-2** swab kits are intended for cleaning the working channels, valves and valve cages of flexible endoscopes.

They are suitable for use with most major brands of gastroscopes, colonoscopes, duodenoscopes, enteroscopes, bronchoscopes, cystoscopes, ureteroscopes, nephroscopes, hysteroscopes and naso-laryngoscopes.

KIT COMPONENTS *(SEE THE TABLE FOR THE DETAILED COMPONENTS BY REFERENCE)*

- **CLEANBRUSH® ENDO**: cleaning brush for working channel
- **CLEANBRUSH® PISTON**: cleaning brush for valves, valve cages
- **CLEANBRUSH® EXO-B**: brush for external cleaning
- **CLEANBRUSH® EXO-S**: sponge for external cleaning

ABOUT THIS MEDICAL DEVICE

1 st release date	1999
Single use	Yes
Class	I
Sterile	No (ethylene oxide decontamination)
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 25 or 50 units depending on product reference
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

PRINCE MEDICAL
 DEVICES

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

CLEANBRUSH® KIT-2

61

REFERENCES FOR CLEANING BRUSHES

AVAILABLE REFERENCE(S)

Kit reference	For working channel	For valves	Packaging
HN1850110	HN1050110 (with 1 brush Ø 1.8 mm L100 cm)	HN1061110 (with 2 brushes Ø11 & 5mm, L15.5 cm)	50
HN1853010	HN1053010 (with 1 brush Ø 3 mm L100 cm)		
HN1853018	HN1053018 (with 1 brush Ø 3 mm L180 cm)		
HN1855010	HN1055010 (with 1 brush Ø 5 mm L100 cm)		
HN1855018	HN1055018 (with 1 brush Ø 5 mm L180 cm)		
HN1855020	HN1055020 (with 1 brush Ø 5 mm L200 cm)		
HN1855023	HN1055023 (with 1 brush Ø 5 mm L230 cm)		
HN1856523	HN1056523 (with 1 brush Ø 6.5 mm L230 cm)		
HN1860110	HN1060110 (with 2 brushes Ø 1.8 mm L100 cm)		
HN1863010	HN1063010 (with 2 brushes Ø 3 mm L100 cm)		
HN1863018	HN1063018 (with 2 brushes Ø 3 mm L180 cm)		
HN1865010	HN1065010 (with 2 brushes Ø 5 mm L100 cm)		
HN1865018	HN1065018 (with 2 brushes Ø 5 mm L180 cm)		
HN1865020	HN1065020 (with 2 brushes Ø 5 mm L200 cm)		
HN1865023	HN1065023 (with 2 brushes Ø 5 mm L230 cm)		
HN1866523	HN1066523 (with 2 brushes Ø 6.5 mm L230 cm)		

KIT COMPONENT REFERENCES

Kit reference	Brush for external cleaning	Sponge for external cleaning	Packaging
HN2750000	HN2580000 (U-shaped with soft hair, useful length 2 x 8 cm)	HN1790000	25

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

* Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

PRINCE MEDICAL
 DEVICES