


CLEANBRUSH® KIT-3

62

3-PART CLEANING SWAB KIT



The **CLEANBRUSH® KIT-3** swab kits are intended for cleaning the working channels, valves and valve cages of flexible endoscopes.

They are suitable for use with the main gastroscopes, colonoscopes, duodenoscopes, enteroscopes, bronchoscopes, cystoscopes, ureteroscopes, nephroscopes, hysteroscopes and naso-laryngoscopes on the market.

KIT COMPONENTS (SEE THE TABLE FOR THE DETAILED COMPONENTS BY REFERENCE)

- **CLEANBRUSH® ENDO**: cleaning brush for working channel
- **CLEANBRUSH® ENDO** (AN1055002, AN1056523): cleaning brush for aspiration channel
- **CLEANBRUSH® ENDO** (AN1050123): cleaning brush for air/water channel, and/or Water Jet and/or balloon inflation channel
- **CLEANBRUSH® PISTON**: cleaning brush for valves, valve cages
- **CLEANBRUSH® EXO-B**: brush for external cleaning

ABOUT THIS MEDICAL DEVICE

1 st release date	1999
Single use	Yes
Class	I
Sterile	No (ethylene oxide decontamination)
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 25 or 50 units depending on product reference
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

PRINCE MEDICAL
 DEVICES

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

CLEANBRUSH® KIT-3

AVAILABLE REFERENCE(S)

REFERENCES FOR CLEANING BRUSHES

OTHER KIT COMPONENTS

63

Kit reference	For working channel	For aspiration channel	For valves	For air/water channel, water jet, balloon inflation channel	External cleaning brush	Packaging
HN1883018	HN1053018 (with 1 brush Ø 3 mm L180 cm)				HN1050123 (with 1 brush Ø 1.8 mm L230 cm)	50
HN1885023	HN1055023 (with 1 brush Ø 5 mm L230 cm)					
HN1886523	HN1056523 (with 1 brush Ø 6.5 mm L230 cm)					
HN1893018	HN1063018 (with 2 brushes Ø 3 mm L180 cm)					
HN1895023	HN1065023 (with 2 brushes Ø 5 mm L230 cm)					
HN1896523	HN1066523 (with 2 brushes Ø 6.5 mm L230 cm)					
HN2770110	HN1050110 (with 1 brush Ø 1.8 mm L100 cm)				HN1061110 (with 2 brushes Ø 11 & 5 mm, L15.5 cm)	50
HN2780110	HN1060110 (with 2 brushes Ø 1.8 mm L100 cm)					
HN2183010	HN1053010 (with 1 brush Ø 3 mm L100 cm)	HN1055002 (with 1 brush Ø 5 mm, L20 cm)				
HN2183018	HN1053018 (with 1 brush Ø 3 mm L180 cm)					
HN2193010	HN1063010 (with 2 brushes Ø 3 mm L100 cm)					
HN2193018	HN1063018 (with 2 brushes Ø 3 mm L180 cm)					
HN2953018	HN1053018 (with 1 brush Ø 3 mm L180 cm)	HN1055023 (with 1 brush Ø 5 mm, L230 cm)				
HN2963018	HN1063018 (with 2 brushes Ø 3 mm L180 cm)					
HN4475018	HN1055018 (with 1 brush Ø 5 mm L180 cm)					
HN4475023	HN1055023 (with 1 brush Ø 5 mm L230 cm)	HN1056523 (with 1 brush Ø 6.5 mm, L230 cm)				
HN4495018	HN1065018 (with 2 brushes Ø 5 mm L180 cm)				HN2580000 (U-shaped with soft hair, useful length 2 x 8 cm)	25
HN4495023	HN1065023 (with 2 brushes Ø 5 mm L230 cm)					
HN2520110	HN1050110 (with 1 brush Ø 1.8 mm L100 cm)					
HN2523007	HN1053007 (with 1 brush Ø 3 mm L75 cm)					
HN2523010	HN1053010 (with 1 brush Ø 3 mm L100 cm)					
HN2523018	HN1053018 (with 1 brush Ø 3 mm L180 cm)					
HN2525010	HN1055010 (with 1 brush Ø 5 mm L100 cm)					
HN2525023	HN1055023 (with 1 brush Ø 5 mm L230 cm)					
HN2525030	HN1055030 (with 1 brush Ø 5 mm L300 cm)					
HN2530110	HN1060110 (with 2 brushes Ø 1.8 mm L100 cm)					
HN2533010	HN1063010 (with 2 brushes Ø 3 mm L100 cm)					
HN2533018	HN1063018 (with 2 brushes Ø 3 mm L180 cm)					
HN2535010	HN1065010 (with 2 brushes Ø 5 mm L100 cm)					
HN2535023	HN1065023 (with 2 brushes Ø 5 mm L230 cm)					
HN2535030	HN1065030 (with 2 brushes Ø 5 mm L300 cm)					

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.