


# CLEANBRUSH® PISTON

## CLEANING BRUSH FOR ENDOSCOPE VALVES

57



The **CLEANBRUSH® PISTON** cleaning brushes are intended for cleaning the valve and valve cages of flexible endoscopes.

They are suitable for use with the main gastroscopes, colonoscopes, duodenoscopes, enteroscopes, bronchoscopes, cystoscopes, ureteroscopes, nephroscopes, hysteroscopes and naso-laryngoscopes on the market.

### AVAILABLE REFERENCE(S)

| Reference | Brush diameter (mm) | Number of brushes | Type of cleaning           | Length (cm) |
|-----------|---------------------|-------------------|----------------------------|-------------|
| HN1061110 | 11 & 5              | 2                 | For valves and valve cages | 15.5        |



### ABOUT THIS MEDICAL DEVICE

|                              |                                             |
|------------------------------|---------------------------------------------|
| 1 <sup>st</sup> release date | 1999                                        |
| Single use                   | Yes                                         |
| Class                        | I                                           |
| Sterile                      | No                                          |
| Contains latex               | No                                          |
| Lifetime                     | 5 years                                     |
| Certification                | CE                                          |
| Packaging                    | Box of 50 units in non-sterile wrapping     |
| Storage                      | In its original box in a cool and dry place |

PRINCE MEDICAL SAS  
 ZA La Sente Du Moulin  
 64 rue Gutenberg - 60530 Ercuis - France  
[www.prince-medical.com](http://www.prince-medical.com)  
 Tel: +33 (0)3 44 26 54 76  
[prince-medical@omerin.com](mailto:prince-medical@omerin.com)

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

PRINCE MEDICAL  
 DEVICES