

ENDOLINE® 3V-TRACT

14

3-LUMENS BALLOON CATHETER



TRANSPARENT
TUBE
AS4560720

The **ENDOLINE® 3V-TRACT** 3-lumens balloon catheter enables the removal of gallstones and injection of contrast material or any other medicated substance into the bile ducts. This medical device is intended mainly for persons with hepatic colic following the appearance of gallstones.

ADVANTAGES

- Product compatible for use with the main endoscopes on the market.
- Radiopaque.
- Standard connector for connecting a syringe for contrast medium injection.
- Stainless steel ring and stainless steel ligatures on the balloon enhance radiopacity of the distal part for better balloon identification.
- The guide path and the contrast path are different.

AVAILABLE BALLOONS



LATEX BALLOON CATHETER

PRECAUTIONS FOR USE
DO NOT USE IN PATIENTS
ALLERGIC TO LATEX



SILICONE BALLOON CATHETER

for patients allergic
to latex.



ABOUT THIS MEDICAL DEVICE

1 st release date	1996
Single use	Yes
Class	Ila
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Latex present	Yes except reference AS3860720
Lifetime	2 years except reference AS3860720 5 years
Certification	CE
Packaging	Box of 5 units in individual sterile wrapping
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS
ZA La Sente Du Moulin
64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

PRINCE MEDICAL
DEVICES

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

ENDOLINE® 3V-TRACT

15

LATEX
BALLOON

Reference	Catheter diameter (mm/Fr)	Inflated balloon diameter (mm)	Balloon inflation max gas (mL)	Balloon injection max liquid (mL)	Wire-guide recommended (inch)	Syringe (cc)	Useful length (cm)	Option
AS1420720	2.35 / 7	16	4.5	2.5	.035	5	200	-
AS4320720	2.35 / 7	16	4.5	2.5	.035	5	200	Transparent
AS4320725	2.35 / 7	16	4.5	2.5	.035	5	200	Transparent
AS4560720	2.35 / 7 1.70 / 5	16	4.5	2.5	.035	5	200	Transparent

SILICONE
BALLOON

AS3860720	2.35 / 7	16	4.5	2.5	.035	5	200	-
-----------	----------	----	-----	-----	------	---	-----	---

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

PRINCE MEDICAL
 DEVICES