26

# **APPLISIL**®SA1A

# STRAIGHT PLUNGER APPLICATOR FOR CREAMS AND GELS



The APPLISIL\* straight plunger applicator is used to apply creams and gels vaginally. It is printable and meets the requirements for legible graduation and regulatory identification of medical dosing devices.

# **CHARACTERISTICS**

Barrel	PE (pharmacopoeia grade) adaptable to M9 tube
Plunger	PS, straight, printable
Colour (standard)	White plunger

#### **AVAILABLE ACCESSORIES**

Protective cap	PE
Use	For straight plunger applicator
Benefit(s)	Improved hygiene

Accessory	Item code	Overall Length (mm)	Diameter (mm)
Сар	1ACH1125	26	11.4

# STANDARD REFERENCE(S)

	٨٩٨١	ICII® CA1A ADDIICA	ATOP				
APPLISIL® SA1A APPLICATOR							
Item code	Volume (ml)	Graduation (ml)	Length (mm)	Tip diameter (mm)			
1ASV1119	0.5	0.5	122	8			
1ASV1128	1	0.5 > 1	122	8			
1ASV1122	2	0.5 > 2	122	8			
1ASV1134	5	1 > 5	122	8			

#### **Union Plastic**

ZA Robert BP 27 43140 Saint-Didier-en-Velay / France www.union-plastic.com Tel.: +33 **(0)4 71 61 13 09** union-plastic@omerin.com

# **OPTIONS**

 Specific colour as per Pantone
Customised printing (logo, graduation, colour)

· Packaging: bulk, flow pack



The information given in this technical data sheet is indicative and subject to change without notice. Union Plastic can under no circumstances be held responsible for any incidents arising from inappropriate use, failure to comply with the rules of the trade and standards in force, or failure to comply with the instructions for use supplied. Can be used without special training or skills (adults in paediatric cases). Non-contractual illustrations and photos. Reproduction prohibited without prior consent from Union Plastic.