

**APPLISIL<sup>®</sup> SAZA****CROSS-SHAPED PLUNGER APPLICATOR FOR CREAMS AND GELS**

27

**VAGINAL ROUTE**

The APPLISIL<sup>®</sup> cross-shaped plunger applicator is used to apply creams and gels vaginally. It offers the highest level of competitiveness without compromising on essential functional features.

**PRODUCT INFORMATION**

Barrel	PE
Plunger	PE
Use	Adaptable to M9 tube
Colour (standard)	White plunger

**STANDARD REFERENCE(S)**

APPLISIL <sup>®</sup> SAZA APPLICATOR			
Item code	Volume (ml)	Length (mm)	Tip diameter (mm)
1ASV0015	1.5	127	8
1ASV0014	2	127	8
1ASV1117	2.5	127	8
1ASV1118	5	127	8
1ASV0025	5	126	9.5

**OPTIONS**

- Specific colour as per Pantone
- Customised printing (logo, graduation, colour)
- Packaging: bulk, flow pack

**Union Plastic**  
 ZA Robert BP 27  
 43140 Saint-Didier-en-Velay / France  
[www.union-plastic.com](http://www.union-plastic.com)  
 Tel.: +33 (0)4 71 61 13 09  
[union-plastic@omerin.com](mailto:union-plastic@omerin.com)



The information given in this technical data sheet is indicative and subject to change without notice. Union Plastic can under no circumstances be held responsible for any incidents arising from inappropriate use, failure to comply with the rules of the trade and standards in force, or failure to comply with the instructions for use supplied. Can be used without special training or skills (adults in paediatric cases). Non-contractual illustrations and photos. Reproduction prohibited without prior consent from Union Plastic.