

APPLISIL[®] SA5A**CROSS-SHAPED PLUNGER APPLICATOR FOR TABLETS**

29

VAGINAL ROUTE

The APPLISIL[®] cross-shaped Plunger applicator is used to apply tablets and ovules vaginally. Different versions can be designed according to tablet size.

PRODUCT INFORMATION

Barrel	PP (pharmacopoeia grade)
Plunger	PE (pharmacopoeia grade). Cross-shaped
Colour (standard)	White plunger
TABLET SIZE:	+/- L 25 x W 10 x H 7 mm

STANDARD REFERENCE(S)

APPLISIL [®] SA5A APPLICATOR		
Item code	Length (mm)	Tip diameter (mm)
1AAV1137*	135	14

*Product in development (for more information, please contact us)

**OPTIONS**

- Specific colour as per Pantone
 - Customised printing (logo, graduation, colour)
- Packaging: bulk, flow pack

Union Plastic
 ZA Robert BP 27
 43140 Saint-Didier-en-Velay / France
www.union-plastic.com
 Tel.: +33 (0)4 71 61 13 09
union-plastic@omerin.com



The information given in this technical data sheet is indicative and subject to change without notice. Union Plastic can under no circumstances be held responsible for any incidents arising from inappropriate use, failure to comply with the rules of the trade and standards in force, or failure to comply with the instructions for use supplied. Can be used without special training or skills (adults in paediatric cases). Non-contractual illustrations and photos. Reproduction prohibited without prior consent from Union Plastic.