

APPLISIL[®] SA1A26 **STRAIGHT PLUNGER APPLICATOR FOR CREAMS AND GELS****VAGINAL ROUTE**

The APPLISIL[®] straight plunger applicator is used to apply creams and gels vaginally. It is printable and meets the requirements for legible graduation and regulatory identification of medical dosing devices.

CHARACTERISTICS

Barrel	PE (pharmacopoeia grade) adaptable to M9 tube
Plunger	PS, straight, printable
Use	Adaptable to M9 tube
Colour (standard)	White plunger

AVAILABLE ACCESSORIES

Protective cap	PE
Use	For straight plunger applicator
Benefit(s)	Improved hygiene



PROTECTIVE CAP

Accessory	Item code	Overall Length (mm)	Diameter (mm)
Cap	1ACH1125	26	11.4

STANDARD REFERENCE(S)

APPLISIL [®] SA1A APPLICATOR				
Item code	Volume (ml)	Graduation (ml)	Length (mm)	Tip diameter (mm)
1ASV1119	0.5	0.5	122	8
1ASV1128	1	0.5 > 1	122	8
1ASV1122	2	0.5 > 2	122	8
1ASV1134	5	0.5 > 5	122	8

OPTIONS

- Specific colour as per Pantone
- Customised printing (logo, graduation, colour)
- Packaging: bulk, flow pack

Union Plastic
 ZA Robert BP 27
 43140 Saint-Didier-en-Velay / France
www.union-plastic.com
 Tel.: +33 (0)4 71 61 13 09
union-plastic@omerin.com



The information given in this technical data sheet is indicative and subject to change without notice. Union Plastic can under no circumstances be held responsible for any incidents arising from inappropriate use, failure to comply with the rules of the trade and standards in force, or failure to comply with the instructions for use supplied. Can be used without special training or skills. Non-contractual illustrations and photos. Reproduction prohibited without prior consent from Union Plastic.