

VAGINAL ROUTE

The APPLISIL® cross-shaped plunger applicator is used to apply creams and gels vaginally. It offers the highest level of competitiveness without compromising on essential functional features.

PRODUCT INFORMATION

Barrel	PE (pharmacopoeia grade)
Plunger	PE (pharmacopoeia grade)
Use	Adaptable to M9 tube
Colour (standard)	White plunger

STANDARD REFERENCE(S)

APPLISIL® SA2A APPLICATOR				
ltem code	Volume (ml)	Length (mm)	Tip diameter (mm)	
1ASV0015	1.5	127	8	
1ASV0014	2	127	8	
1ASV1117	2.5	127	8	
1ASV1118	5	127	8	
1ASV0025	5	126	9.5	

Union Plastic

ZA Robert BP 27 43140 Saint-Didier-en-Velay / France www.union-plastic.com Tel.: +33 **(0)4 71 61 13 09** union-plastic@omerin.com

OPTIONS

 Specific colour as per Pantone
Customised printing (logo, graduation, colour)
Packaging: bulk, flow pack



The information given in this technical data sheet is indicative and subject to change without notice. Union Plastic can under no circumstances be held responsible for any incidents arising from inappropriate use, failure to comply with the rules of the trade and standards in force, or failure to comply with the instructions for use supplied. Can be used without special training or skills. Non-contractual illustrations and photos. Reproduction prohibited without prior consent from Union Plastic.