## **FOR 1-HOLE PIPETTE**



**5AA** adapter
For pipette with rounded tip
and PP28 bottleneck



**5CA** adapter For connect pipette and PP28 bottleneck

# FOR ATRAUMATIC 3-HOLES PIPETTE



**5TA** adapter For atraumatic pipette

### FOR DOSING SYRINGE



**6CA** adapter with fins



**6CA** adapter with short tube



**6CA** adapter with long tube

#### Union Plastic

ZA Robert BP 27 43140 Saint-Didier-en-Velay / France www.union-plastic.com Tel.: +33 **(0)4 71 61 13 09** union-plastic@omerin.com



# DOS'UP®5AA, 5CA, 5TA, 6CA

*ADAPTERS* 

17

# ORAL ROUTE

**DOS'UP®** adapters are used with DOS'UP® pipettes and dosing syringes. They serve to optimise dosing with improved hygiene.

### PRODUCT INFORMATION

Adapters

Material	PE (food and/or pharmacopoeia grade)
	For PP28 bottleneck, glass and plastic (see reference table).
	> Models for 1-hole pipette with Rounded or Connect tip
Use	> Model for 3-holes pipette
	> Models for dosing syringes

## STANDARD REFERENCE(S)

DOS'UP® PIPETTE ADAPTER							
Item code	Reference	Bottleneck diameter	Vial type	Adapter shape	Pipette type		
1DAD0202	5AA	PP28	Glass	Fins	1 hole		
1DAD0203	5AA	PP28	Plastic	Fins	1 hole		
20034500*	5CA	PP28	Glass and plastic	Short tube	Connect		
1DAD1110	5TA	PP28	Glass and plastic	Fins	3 holes		
1DAD1111	5TA	PP28	Glass and plastic	Fins	3 holes		
1DAD1112	5TA	PP28	Glass and plastic	Fins	3 holes		
1DAD1113	5TA	PP28	Glass and plastic	Fins	3 holes		
1DAD1114	5TA	PP28	Glass and plastic	Fins	3 holes		

DOS'UP® DOSING SYRINGE ADAPTER								
Item code	Reference	Bottleneck diameter	Vial type	Adapter shape	Dosing syringe			
210004500*	6BA	GL18	Glass	Short tube	0.2 and 1 ml			
1DAD0241	6CA	PP28	Glass and plastic	Fins	5.6, 10 and 12 ml			
1DAD0239	6CA	PP28	Glass	Short tube	5.6, 10 and 12 ml			
1DAD1104	6CA	PP28	Glass	Long tube	5.6, 10 and 12 ml			

<sup>\*</sup>Product in development (for more information, please contact us)

The information given in this technical data sheet is indicative and subject to change without notice. Union Plastic can under no circumstances be held responsible for any incidents arising from inappropriate use, failure to comply with the rules of the trade and standards in force, or failure to comply with the instructions for use supplied. Can be used without special training or skills (adults in pediatric cases). Non-contractual illustrations and photos. Reproduction prohibited without prior consent from Union Plastic.