

# MEDIC'UP<sup>®</sup> VA2A

## VIAL TO VIAL CONNECTING DEVICE

49

### PARENTERAL ROUTE

The MEDIC'UP<sup>®</sup> vial to vial connecting device makes it quick and easy to connect two Vial Adapters.

#### CHARACTERISTICS

Spike	Made of plastic, allowing the rubber stopper penetration without deformation or chipping
3 Holes	> Facilitates the transfer of mixtures between two bottles > Optimises rinsing, liquid transfer and aspiration of the reconstituted medicine
Luer-Lock	ISO 80369-compatible male and female connectors > Guide system allowing perfect central penetration of the tip at the level of the vial's rubber stopper
Benefit(s)	> Compatible with all Ø20 and Ø13 bottlenecks > Permanent clip function on the bottleneck > Moulded in Polycarbonate > Sterilised

#### STANDARD REFERENCE(S)

MEDIC'UP <sup>®</sup> VIAL TO VIAL		
Item code*	Diameter of male Luer vial (mm)	Diameter of female Luer vial (mm)
23602000*	20	20/13

\*Product in development (for more information, please contact us)



**Union Plastic**  
 ZA Robert BP 27  
 43140 Saint-Didier-en-Velay / France  
[www.union-plastic.com](http://www.union-plastic.com)  
 Tel.: +33 (0)4 71 61 13 09  
[union-plastic@omerin.com](mailto:union-plastic@omerin.com)



The information given in this technical data sheet is indicative and subject to change without notice. Union Plastic can under no circumstances be held responsible for any incidents arising from inappropriate use, failure to comply with the rules of the trade and standards in force, or failure to comply with the instructions for use supplied. Medical devices marketed by Union Plastic must only be used by a duly qualified practitioner or authorised person.  
 ® Registered trademark of Union Plastic. Non-contractual illustrations and photos. Reproduction prohibited without prior consent from Union Plastic.