

# VETILINE® IM1A

## INTRAMAMMARY SYRINGE

27

### INTRAMAMMARY ROUTE

The VETILINE® intramammary syringe is ideal for injecting drugs into cows' udders thanks to its extremely fine needle. The snap-off cap allows easy use.

#### PRODUCT INFORMATION

- Syringe
- > Food and/or pharmacopoeia grade
  - > Barrel with fitted cap
  - > Regular or snap-off cap (see table below)
  - > Separate straight plunger

Colour (standard)	White (specific colour as per Pantone: please contact us)
Benefits	> Easy use (references with snap-off cap) > Dosing accuracy (Ø 1.5 mm fine cannula)
Option	Optional digital labelling

#### STANDARD REFERENCE(S)

##### VETILINE® IM1A SYRINGE WITH REGULAR CAP

Item code Barrel + cap	Volume (ml)	Cannula Diameter (mm)	Barrel Diameter (mm)	Length Barrel + cap (mm)	Item code Plunger	Plunger Type	Plunger Length (mm)
1TSE0654	10	1.5	20.5	101	1TPI0640	Straight	84
1TSE0655	20	1.5	22	113	1TPI0757	Straight	89
1TSE0684	30	1.5	26	134	1TPI0750	Straight	113

##### VETILINE® IM1A SYRINGE WITH SNAP-OFF CAP

1TSE0661	5	1.5	14	75	1TPI0761	Straight	43
1TSE1127	7	1.5	13	90	1TPI0639	Straight	70
1TSE1300	12	1.5	19	103	1TPI0642	Straight	88

SNAP-OFF CAP



NF EN ISO 15378: 2018



SMOOTH PLUNGER

BARREL

SNAP-OFF CAP



**Union Plastic**  
 ZA Robert BP 27  
 43140 Saint-Didier-en-Velay-France  
[www.union-plastic.com](http://www.union-plastic.com)  
 Tel.: +33 (0)4 71 61 13 09  
[union-plastic@omerin.com](mailto:union-plastic@omerin.com)



The information given in this technical data sheet is indicative and subject to change without notice. Union Plastic can under no circumstances be held responsible for any incidents arising from inappropriate use, failure to comply with the rules of the trade and standards in force, or failure to comply with the instructions for use supplied. Can be used without special training or skills. Non-contractual illustrations and photos. Reproduction prohibited without prior consent from Union Plastic.