


PM-CARE® ENDOBIOPS

HISTOLOGICAL SAMPLING CATHETER



The PM-CARE® ENDOBIOPS histological sampling catheter is used to take histological samples from the endometrium.

ADVANTAGES

- Centimetric graduations for better control of catheter insertion.
- Atraumatic tip enabling smooth introduction into the cervical canal.
- Rigid and notched grip tip for a firm grip.
- Efficient aspiration provided by negative pressure.

AVAILABLE REFERENCE(S)

Reference	Model	Diameter (mm / Fr)	Length (cm)
PS1200923	STANDARD	3 / 9	23
PS1210626	MINI	2 / 6	27
PS1220627	CURVED	2 / 6	27
PS3490923	PLUS	3 / 9	23

ABOUT THIS MEDICAL DEVICE

1 st release date	2001
Single use	Yes
Class	Is
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 25 units in individual sterile wrapping
Storage	In its original box in a cool and dry place



AVAILABLE MODEL(S)



STANDARD



MINI

- For patients with a narrow cervix or cervical canal
- For post-menopausal patients



CURVED

- For swabbing from the lateral side of the endometrium
- For patients with a non-conventional shaped uterus



PLUS

- To improve swabbing quality
- To take two separate swabbing