ENDOLINE® COLD-LOOP

32 COLD SNAR



The ENDOLINE® COLD-LOOP cold snare is intended for endoscopic resection of small sessile/flat (≤ 9 mm) polyps present in the gastrointestinal tract via an upper or lower endoscopy without use of diathermic energy.

ADVANTAGES

- Three fingers grip for easy manipulation.
- Radiopaque distal part for perfect visibility.
- Product compatible for use with the main endoscopes on the market.

AVAILABLE REFERENCE(S)

| Reference | Sheath diameter (mm) | Sheath length (cm) | Opening (mm) |
|-----------|----------------------|-----------------------|-----------------|
| AS4372523 | 2.35 | 230 | 10 x 25 |

ABOUT THIS MEDICAL DEVICE

| 1st release date | 2016 |
|----------------------|-----------------------------------------------|
| Single use | Yes |
| Class | lla |
| Sterile | Yes |
| Sterilisation method | Ethylene oxide sterilisation |
| Contains latex | No |
| Lifetime | 5 years |
| Certification | CE |
| Packaging | Box of 5 units in individual sterile wrapping |
| Storage | In its original box in a cool and dry place |
| | |

PRINCE MEDICAL SAS
ZA La Sente Du Moulin
64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com



The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.