

APPLISIL[®] TU3A

30

SINGLE-DOSE TUBE

VAGINAL, CUTANEOUS OR RECTAL ROUTE



The APPLISIL[®] single-dose tube is used to apply creams and gels vaginally or rectally. Practical and hygienic, this product also has a snap-off tip for easy opening before use. Easily filled, its thin wall allows closure by hot air welding.

CHARACTERISTICS

Material	PE (pharmacopoeia grade)
Tip	Snap-off
Colour (standard)	White
Printable	1 colour

STANDARD REFERENCE(S)

APPLISIL [®] TU3A SINGLE-DOSE TUBE		
Item code	Volume (ml)	Length (mm)
17031400*	1	60
1ATU0002	5	110

*Product in development (for more information, please contact us)

OPTIONS

- Specific colour as per Pantone
- Customised printing (logo, graduation, colour)
- Packaging: bulk

Union Plastic
 ZA Robert BP 27
 43140 Saint-Didier-en-Velay / France
www.union-plastic.com
 Tel.: +33 (0)4 71 61 13 09
union-plastic@omerin.com

Union plastic

The information given in this technical data sheet is indicative and subject to change without notice. Union Plastic can under no circumstances be held responsible for any incidents arising from inappropriate use, failure to comply with the rules of the trade and standards in force, or failure to comply with the instructions for use supplied. Can be used without special training or skills (adults in paediatric cases). Non-contractual illustrations and photos. Reproduction prohibited without prior consent from Union Plastic.