

# MEDIC'UP<sup>®</sup> VA3A

50

TRANSFER DEVICE VIAL TO IV BAG

PARENTERAL ROUTE



The MEDIC'UP<sup>®</sup> transfer device vial to IV Bag is a needleless transfer device for connecting a vial of medication to an IV bag with a rubber inlet.

## CHARACTERISTICS

Double spike	<ul style="list-style-type: none"> <li>&gt; Plastic for chip-free, distortion-free penetration of the rubber stopper on the vial and infusion bag</li> </ul>
Benefit(s)	<ul style="list-style-type: none"> <li>&gt; Guide system allowing perfect central penetration of the tip at the level of the vial's rubber stopper</li> <li>&gt; Facilitates mixture transfer and return</li> <li>&gt; Optimises liquid transfer, rinsing and aspiration of the reconstituted medicine</li> <li>&gt; Compatible with all Ø20 bottlenecks</li> <li>&gt; Permanent clip function on the bottleneck to ensure single use of the transfer and a good hold during use</li> <li>&gt; Sterilised by Gamma irradiation</li> </ul>

## STANDARD REFERENCE(S)

MEDIC'UP <sup>®</sup> VIAL TO IV BAG		
Item code	Vial diameter (mm)	Height (mm)
23603000*	20	34

\*Product in development (for more information, please contact us)

**Union Plastic**

ZA Robert BP 27

43140 Saint-Didier-en-Velay / France

[www.union-plastic.com](http://www.union-plastic.com)

Tel.: +33 (0)4 71 61 13 09

[union-plastic@omerin.com](mailto:union-plastic@omerin.com)



The information given in this technical data sheet is indicative and subject to change without notice. Union Plastic can under no circumstances be held responsible for any incidents arising from inappropriate use, failure to comply with the rules of the trade and standards in force, or failure to comply with the instructions for use supplied. Medical devices marketed by Union Plastic must only be used by a duly qualified practitioner or authorised person.  
 ® Registered trademark of Union Plastic. Non-contractual illustrations and photos. Reproduction prohibited without prior consent from Union Plastic.